


Evaluating the Infectivity and Prognosis of a Virion Sample using Plasmonic Nanoparticles and MATLAB Programming

Bryan Hong




Status Quo

2 Standards of testing:



Nucleic Acid Test



Antigen Test

- **Nucleic Acid Test advantages:** Accurate, sensitive
- **Nucleic Acid Test disadvantages:** Resources, slow


- **Antigen Test advantages:** Rapid, cheap
- **Antigen Test disadvantages:** Not as sensitive, false-positives

Status Quo

Viral Load Quantification:

- Method of counting virions present within a nasal swab
- Predicts the prognosis, infectibility, and severity of the viral sample
- Conducted using **Cycle Threshold Values** via PCR tests

Figure created by Tyson Holmes, 2021



- Counts the number of viral amplification steps before the virus is detected by the PCR test

- **Contains many issues**

Engineering Goals

I. Develop a rapid, sensitive, and accurate RSV diagnostic test that combines the advantages of both the antigen test and the nucleic acid test

II. Create a viral load quantification system that counts and displays the individual number of virions

Principals and Procedures

I.

RSV AuNP antibody aggregation






II.

Nanobubble detection

III.

Virion counting mechanism

Principal I- AuNP Antibody Conjugation




Gold Nanoparticle
(AuNP)

DTSSP
Crosslinker

Palivizumab
Monoclonal
Antibody


Principal I- RSV AuNP Aggregation

Figure
created by
Joseph
Horus, 2020




Palivizumab Antibodies
binding to the F-
glycoproteins on the RSV
Virus

One virion has multiple
AuNPs attached to it




Unaggregated
AuNP-Antibodies



RSV-aggregated
AuNP-Antibodies


Principal II- Nanobubble Detection



Specific light wavelength will cause electrons on NP to oscillate → **Plasmon Resonance**

Plasmon Resonance for 15 nm AuNP = 532 nm (green light)

Principal II- Nanobubble Detection




When the electrons oscillate on the AuNP, the AuNP heats up fluid around it → **nanobubble**

Using a red HeNe laser to detect **nanobubble**

Figure created by Jackson Martin, 2019


Procedure II- Nanobubble Detection



Step 1:

- Green picosecond pulsating laser → energizes the AuNPs → generates **nanobubbles**
- Red HeNe laser → detects the **nanobubbles**

Procedure II- Nanobubble Detection



Step 2:

- Red HeNe laser diffracted by the **nanobubble**
- Laser energy decreases


Procedure II- Nanobubble Detection

Step 3:


- Displayed onto the oscilloscope

X-axis: Time (ns)


Y-axis: Energy of lasers (mV)



Procedure II- Nanobubble Detection



No **nanobubble** detected



Nanobubble detected

Nanobubble
detected →
laser energy
decreases

Principal III- Virion Counting System

I.


Larger AuNP
aggregates → larger
nanobubbles

II.

Larger **nanobubbles**
→ more laser
diffraction


Procedure III- Virion Counting System

Off
signal



Unaggregated
Nanobubble detected

On
signal



RSV-aggregated
Nanobubble detected

**Unaggregated
nanoparticles**



**Smaller
nanobubbles**



**Diffracting less
laser energy**



**Cause smaller dip
on oscilloscope**

**Aggregated
nanoparticles**



**Larger
nanobubbles**



**Diffracting more
laser energy**



**Cause larger dip
on oscilloscope**

Procedure III- Virion Counting System

Step 1:

- Feed the **unaggregated AuNPs** through the microcapillary and record amplitude and AUC of every signal
- Establishes the **threshold**

$$\text{Threshold} = \mu + 5\sigma$$


Amplitude threshold
(T_{AMP})

AUC threshold (T_{AUC})


Procedure III- Virion Counting System

Step 2:

- Feed the **RSV-aggregated AuNPs** into the microcapillary
- Record amplitude and AUC of every signal
- $> T_{AMP}$ or $T_{AUC} \rightarrow$ counted RSV virion




Results



$$T_{AMP} = 225.65$$

$$T_{AUC} = 137.13$$



Conclusions

I. Develop a rapid, sensitive, and accurate RSV diagnostic test that combines the advantages of both the antigen test and the nucleic acid test

- Generated results → avg 5.2 mins
- Accuracy +/- → 99.99%

II. Create a viral load quantification system that counts and displays the individual number of virions

- Results can be cross-validated using **Poisson Statistics**

Impacts







Figure created by Giovanni Piedimonte, 2014

- Allow scientists to gain more information about spread of viruses
- Gives more information to healthcare providers to personalize treatments
- Faster and more accurate diagnostic results help limit spread and increase prevention



$$\lambda = -\ln(1 - f_{on}) \quad (1)$$

$$V = c/\lambda \quad (2)$$




DTSSP

3,3'-Dithiobis(sulfosuccinimidylpropionate)


MW 608.51

Spacer Arm 12.0 Å

- Disulfide bonds get reduced into 2 **thiol** bonds
- Breaks the DTSSP into two parts
- Bonds to the gold ions using **thiol-gold** bonding



- Sulfo-NHS-ester bonds bind to the amine bonds to form amide bonds



- Binds to the heavy chain of antibody

Error Analysis

- Concentration of RSV abnormally high

- Laboratory mistakes

- MATLAB counting errors



Thank You

**All diagrams
were created by
the researcher
unless specified
otherwise**

**All images were
taken by the
researcher
unless specified
otherwise**